# OpenLAB CDS 2.x

## Thai Manual for 8890 GC

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**Front Inlet : Leak & Restriction Test** ? Close

Testing Primary sample flow path.

Description 1. Verify the inlet can control to a pressure setpoint, validating the basic operation of the inlet and flow path.  
2. Hold the pressure and monitor the flow error between actual and target column flow. If the actual is larger than the target, a leak exists in the flow path. If the actual is less than the target, a restriction exists in the flow path.

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
Close Test will update the tests info in the System Health Report

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The screenshot displays the GC8890 - Acquisition software interface. The top navigation bar includes File, Home, GC Plugins, and Control. The Control panel features buttons for Take, Release, Status, Method, Single Sample, Sequence, Copy, Delete, and Reset. The main area is divided into two sections: Instrument Status and Online Signals.

**Instrument Status:** This section shows the status of the GC8890 instrument and its components. The instrument is currently in an 'Idle' state. The Front Inlet Flow Path is also shown, with parameters for Front Inlet (0.3 psi [28.2 psi], 48.7 °C [50 °C]), Column #1 (40 °C [40 °C], 0 mL/min), and Front FID (50.1 °C [50 °C]). A warning message indicates: '\* One or more heated zones are off. This can result in large differences between actuals and setpoints for flow and pressure.'

**Online Signals:** This section displays a 'Test Plot' graph. The y-axis is labeled 'Signal Selection' and ranges from -1 to 1. The x-axis is labeled 'Time (minutes)' and ranges from -52 to 6. The plot shows a single horizontal line at approximately 0.25 signal level.

This screenshot shows the same GC8890 - Acquisition software interface, but with the 'Method' button in the Control panel highlighted in blue, indicating it is the selected option.

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
## Directory Structure of CDS 2.x Software

- Master methods reside in **Project\Methods** folder.
- Result set methods are copies of master methods and reside in **Project\Results\MyResult.rslt** folder.

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

## Create a new Acquisition method

Method has been modified.  
Save before continuing?

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

Edit Properties of Installed GC Column

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Inventory... Create Layout Update

Find Clear  Only Show Column

Drag a column header here to group by that column

Part Number	M...	Description	Fav...	Len... m	Dia... µm	Film Thi...	Phase Ratio	Min Temp,	Max Temp,	Max Prog	Form Factor	Key...	Comments	Time Stamp
▶ CP7527	Agilent	CP-Carbowax...	<input type="checkbox"/>	50	320	0.2	400	60	80	80	7-inch	CP-C...		01/11/19 ...
⊞ CP7494	Agilent	CP-Chirasil-D...	<input type="checkbox"/>	25	250	0.08	781	-1	200	200	7-inch	CP-C...		01/11/19 ...
⊞ CP7495	Agilent	CP-Chirasil-L ...	<input type="checkbox"/>	25	250	0.12	521	-1	200	200	7-inch	CP-C...		01/11/19 ...
⊞ CP749515	Agilent	CP-Chirasil-L ...	<input type="checkbox"/>	25	250	0.12	521	-1	281	281	5-inch	CP-C...		01/11/19 ...
⊞ CP7502	Agilent	CP-Chirasil-D...	<input type="checkbox"/>	25	250	0.25	250	-1	200	200	7-inch	CP-C...		01/11/19 ...

GC Column Catalog x

Inventory... Create Layout Update

Find Clear  Only Show Column

Drag a column header here to group by that column

Part Number	Ma...	Description	Fav...	Len... m	Dia... µm	Film Thi...	Phase Ratio	Min Temp,	Max Temp,	Max Prog	Form Factor	Keyw...	Comments	Time Stamp
▶ 19091J-413	Agilent	HP-5	<input checked="" type="checkbox"/>	30	320	0.25	320	-60	325	350	7-inch	HP-5		01/11/19...
⊞ 19091J-4...	Agilent	HP-5	<input type="checkbox"/>	30	320	0.25	320	-60	325	350	5-inch	HP-5		01/11/19...
⊞ 19091J-4...	Agilent	HP-5	<input type="checkbox"/>	30	320	0.25	320	-60	325	350	LTM	HP-5...		01/11/19...
⊞ 19091J-4...	Agilent	HP-5	<input type="checkbox"/>	30	320	0.25	320	-1	281	281	Intuvo	Polys...		01/11/19...

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Average Velocity: 85.075 cm/sec

Holdup Time: 0.58772 min

(Initial): 0 min  
He @ 40 °C Oven  
Out: Ambient Pressure  
30 m x 320 µm x 0.25 µm

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Column Outlet Pressure: 0 psi

	Column	Calibration Results	Inlet	Outlet	Heated By
1	Agilent 19091S-433UI: <Not Invented> HP-5ms Ultra Inert -60 °C—325 °C (350 °C): 30 m x 250 µm x 0.25 µm	Uncalibrated	Front Inlet	Front Det...	Oven
2	Agilent 19091J-413: <Not Invented> HP-5 -60 °C—325 °C (350 °C): 30 m x 320 µm x 0.25 µm	Uncalibrated	Back Inlet	Front Det...	Oven
3	No Column Installed	Uncalibrated	Unspecified	Other	Oven
4	No Column Installed	Uncalibrated	Unspecified	Other	Oven
5	No Column Installed	Uncalibrated	Unspecified	Other	Oven
6	No Column Installed	Uncalibrated	Unspecified	Other	Oven

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Acquisition Method – Untitled

The screenshot displays the 'Acquisition Method – Untitled' window. On the left, there are navigation tabs for 'General Properties' and 'Instrument Setup Agilent 7890B'. The main area is divided into three sections: a flow path diagram, a graph, and a configuration panel. The flow path diagram shows two paths: 'Front Inlet Flow Path' and 'Back Inlet Flow Path', each with a 'Front SS Inlet' and 'Column #1' and 'Column #2'. The graph plots 'Oven: °C\*' (red line), 'Column #2 He: psi' (black line), and 'Column #2 He: mL/min' (green line) against 'Run Time, min'. The configuration panel includes 'Injection' settings (Syringe Size: 10 µL, Injection Volume: 1 µL) and 'Washes and Pumps' settings (Preinj, Postinj, Volume (µL) for Solvent A, Solvent B, and Sample Washes, and Sample Pumps).

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Split-Splitless Inlet

Select Liner... Options not been selected.

Parameter	Actual	Setpoint
Heater	50 °C	50 °C
Pressure	.091 psi	29.063 psi
Total Flow	-0.1172 mL/min	85.236 mL/min
Septum Purge Flow	0.3712 mL/min	3 mL/min

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### Acquisition Method – Untitled


The screenshot shows the 'Acquisition Method – Untitled' window for an Agilent 7890B. The 'Signals' tab is active, displaying a 'Signal Event Table' and a 'Dual' configuration table. The 'Dual' table is as follows:

Dual	Signal Source	Data Rate / Min Peak Width	Zero	Save
F	#1: Front Signal (FID)	50 Hz / 0.004 min	Hz ?	<input checked="" type="checkbox"/>
B	#2: None	50 Hz / 0.004 min	Hz ?	<input type="checkbox"/>
B	#3: None	50 Hz / 0.004 min	Hz ?	<input type="checkbox"/>
B	#4: None	50 Hz / 0.004 min	Hz ?	<input type="checkbox"/>


Below the table is a 'Signal Event Table' with columns for 'Signal Source', 'Time, min', and 'Signal Event'. A 'Delete Events' button is also visible.

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


### 1. Click Single Sample

The screenshot shows the GC8890 - Acquisition software interface. The top menu bar includes File, Home, GC Plugins, and Control. The 'Single Sample' button is highlighted in the 'Layouts' section. The 'Single Sample Analysis' window is open, displaying the following fields:

- Sample name: [ ]
- Acq. method: C:\CDSPProject\Training\Method\Training.amx
- Proc. method: [ ]
- Result path: C:\CDSPProject\Training\Result
- Result name: [ ]

The 'Run' button is highlighted in green. The 'Run Queue' window is open, showing a table with columns: State, Type, Result Name, User, Acquisition Met, and Details. The 'Online Signals' window is open, showing a graph for the Agilent 7890B instrument with a 'Front Signal' plot. The y-axis ranges from -1 to 1, and the x-axis ranges from 186 to 244. The signal is currently at 0.

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
1. Click Sequence

2. กด  (New Sequence)

The screenshot shows the GC8890 - Acquisition software interface. The 'Sequence' tab is active, displaying a table for configuring the sequence run. The table has columns for Action, Vial, Sample type, Run type, Level, Acq. method, Proc. method, Inj/Vial, Volume, Injection source, Sample amount, Sample name, and Data file. The table contains 7 rows of data, including 'Cal. Std.' and 'Sample' entries. Below the table, there are sections for 'Sample custom parameters', 'Compound custom parameters', 'Compound amounts', 'Description', and 'Create report'. At the bottom, there are fields for 'Result path' (C:\CDSPProject\Training\Result) and 'Result name', along with a green 'Run' button.

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
### 1. Double Click OpenLab Control Panel

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The screenshot shows the OpenLab Control Panel interface. At the top, there is a 'MANAGEMENT' section with icons for 'Create', 'Edit Project', 'Delete Project', 'Refresh All', 'Start Data Analysis', and 'Create Desktop Shortcut'. Below this is a 'Projects and Groups' section with a 'Projects' sidebar. The sidebar shows a tree view with 'Projects' expanded, containing 'Agilent', 'PT', and 'Training'. The 'Training' project is selected. The main panel displays the 'Training' project details, including 'Name: Training', 'Project folder path: C:\CDSPProjects\Training', and 'Description:'. There are also expandable sections for 'Details' and 'Activity Log (last 7 days)'. At the bottom, there are tabs for 'Instruments', 'Projects', and 'Administration', with 'Projects' currently selected.

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The screenshot shows the 'Classroom - Data Analysis' interface. On the left, the 'Data Selection' panel lists various data files under the 'Classroom' folder. On the right, the 'Injection List' panel displays a table of injection records.

Order No	Inj. #	Sample name	Data file	Description	Sample type
1	1	Calibration 1.1	_001_008-0101.D		Cal. Std.
2	1	Calibration 1.2	_002_008-0201.D		Cal. Std.
3	1	Calibration 2.1	_003_008-0501.D		Cal. Std.
4	1	Calibration 2.2	_004_008-0701.D		Cal. Std.
5	1	Calibration 3.1	_005_008-1001.D		Cal. Std.
6	1	Calibration 3.2	_006_008-1101.D		Cal. Std.
7	1	Sample 1	_007_008-0301.D		Sample
8	1	Sample 2	_008_008-0401.D		Sample
9	1	Sample 3	_009_008-0801.D		Sample
10	1	Sample 4	_010_008-0901.D		Sample

The screenshot shows the 'Classroom - Data Analysis' interface with the 'Data Processing' panel on the left, a 'Chromatograms' plot in the center, and two windows on the right: 'Sample Information' and 'Injection Results'.

**Chromatograms**  
 The plot shows Response [mAU] vs Retention time [min]. Peaks are labeled with their retention times: 1.353, 4.265, 5.093, 6.687, 6.764, and 9.574. The plot title is 'Calibration 1.1 | DAD1 A, Sig=254.4 Ref=360,100 | \_001\_008-0101.D'.

**Sample Information**

Sequence	
Name	Example_D
Description	Self assem
Creation date (yyyy-MM-dd)	
Created by	SYSTEM
Modification date (yyyy-MM-dd)	2019-10-09
Modified by	SYSTEM
Sample	
Name	Calibration
Description	

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



The screenshot shows the Classroom - Data Analysis software interface. The top menu bar includes File, Home, Processing, Audit/E-Sign, and Injection Tree. The main toolbar contains icons for Open Data, Close Data, Open Method, Save Method, Close Method, Reprocess All, Save All Results, Print All, View PDF, and a Windows icon. A red box highlights the Windows menu, which includes Acquisition Setpoints, Injection List, Peak Details, Calibration Curve, Injection Results, Sample Information, Chromatograms, Processing Method, Peak Explorer, and Windows.

The Data Processing window is open, showing a list of data files under 'Example\_Data\_for\_ISTD':

- ISTD-level1 - \_001\_ISTD-LV-000001.D
- ISTD-level2 - \_002\_ISTD-LV-000002.D
- ISTD-level3 - \_003\_ISTD-LV-000003.D
- ISTD-sample - \_004\_ISTD-SAMPI F-001.D

The Processing Method window is open, showing the 'Example\_Method\_for\_ISTD' method. The 'General' section is selected, and the 'Integration Events ChemStation' section is expanded to show 'Standard' and 'Advanced' options. The 'Signals' section is also visible, showing a list of signals and their corresponding times.

The Chromatograms window is open, showing a chromatogram plot with 'Response [mAU]' on the y-axis and 'Retention time [min]' on the x-axis. The plot shows two peaks labeled 'Nifedipin' at 4.190 and 4.892 minutes, and another peak labeled 'Nitrendipin' at 5.684 minutes. The plot title is 'ISTD-level1 | DAD1 A, Sig=235,10 Ref=550,60 | \_001\_ISTD-LV-000001.D'.

The Injection List window is open, showing a table with columns for Sample a, ISTD Amt 1, ISTD Amt 2, ISTD Amt 3, and ISTD Amt 4. The table contains four rows of data, all with a value of 0.50000 in the ISTD Amt 1 column.

Y1	Sample a	ISTD Amt 1	ISTD Amt 2	ISTD Amt 3	ISTD Amt 4
..	0.00000	0.50000			
..	0.00000	0.50000			
..	0.00000	0.50000			
..	0.00000	0.50000			

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1. Start Data Analysis > Data selection > Load data > Data Processing

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The screenshot displays the 'Classroom - Data Analysis' software interface. The top ribbon shows the 'Processing' tab selected (1). The 'New Method' button in the ribbon is highlighted (2). Below the ribbon, the 'Data Processing' pane shows a list of data files, and the 'Chromatograms' pane shows a plot. The 'Create New Processing Method' dialog box is open, showing a list of method configurations. The 'GC/LC Quantitative' option is selected (3). The 'Create method' button is highlighted (4).

**Create New Processing Method**

Select method configuration

- GC/LC Area Percent
- GC/LC Quantitative**
- 3D UV Quantitative
- MS
- LC/MS Sample Purity

GC/LC Quantitative

This method configuration supports:


- Integration & Area%/Height% calculation
- Signal alignment
- Blank subtraction for GC, LC, A/D data
- Compound identification
- Calibration
- System Suitability
- Custom Calculation
- Post Processing Plugins
- Reporting

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The screenshot displays the 'Classroom - Data Analysis' software interface. The top ribbon includes tabs for File, Home, Processing, Audit/E-Sign, and Injection Tree. The 'Processing' tab is active, showing options like 'Reprocess All', 'Save All Results', 'Print All', and 'View PDF'. The 'Chromatograms' and 'Processing Method' buttons are highlighted in the ribbon.

The main window is divided into three panes:

- Data Processing:** Shows a tree view under 'Example\_Data\_for\_ISTD' with four sub-items: 'ISTD-level1 - \_001\_ISTD-LV-000001.D', 'ISTD-level2 - \_002\_ISTD-LV-000002.D', 'ISTD-level3 - \_003\_ISTD-LV-000003.D', and 'ISTD-sample - \_004\_ISTD-SAMPLE-001.D'. Each item has a checkmark icon.
- Processing Method:** Shows the configuration for the 'Police' method. The 'Info' tab is selected, displaying:
  - Location: C:\CDSPr...Police.pmx
  - Description: (Empty text box)
  - Created at (yyyy-MM-dd): 2019-10-10 21:41:37+07:00
  - Modified at (yyyy-MM-dd): 2019-10-10 21:41:37+07:00
- Chromatograms:** Shows a grid for plotting. The y-axis is labeled 'Response' and ranges from 1.5 to 5.0 (scaled by 10<sup>2</sup>). The x-axis is unlabeled. A message box in the center of the grid states 'No injection selected'.

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The screenshot shows the 'Processing Method' window with the 'Standard' integration event selected. The 'Chromatograms' window displays a chromatogram with peaks at 4.190, 4.682, 5.654, and 6.701 minutes.

Use	Time (min)	Event	Value
<input checked="" type="checkbox"/>	0.000	Slope sensitivity	1.00000
<input checked="" type="checkbox"/>	0.000	Peak width	0.02000
<input checked="" type="checkbox"/>	0.000	Area reject	1.00000
<input checked="" type="checkbox"/>	0.000	Height reject	1.70000
<input checked="" type="checkbox"/>	0.000	Shoulders mode	Off
<input checked="" type="checkbox"/>	0.000	Area% reject	0.00000

## Processing Method


The detailed view of the 'Processing Method' window shows the 'Standard' integration event configuration. The 'Chromatograms' window displays a chromatogram with peaks at 4.190, 4.682, 5.654, and 6.701 minutes.

Use	Time (min)	Event	Value
<input checked="" type="checkbox"/>	0.000	Slope sensitivity	1.00000
<input checked="" type="checkbox"/>	0.000	Peak width	0.02000
<input checked="" type="checkbox"/>	0.000	Area reject	1.00000
<input checked="" type="checkbox"/>	0.000	Height reject	1.70000
<input checked="" type="checkbox"/>	0.000	Shoulders mode	Off
<input checked="" type="checkbox"/>	0.000	Area% reject	0.00000

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
Integration Events ChemStation > Standard > Global, right-click next to table > Add integration

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


6. กด Print All > View PDF

## Single Injection Report Agilent Technologies

**Sample name:** ISTD-sample  
**Data file:** \_004\_ISTD-SAMPLE-001.D  
**Instrument:** Instrument 1  
**Inj. volume:** 2.000  
**Acq. method:** Example\_Method\_for\_ISTD.pmx  
**Processing method:** Example\_Method\_for\_ISTD.pmx  
**Manually modified:** None

**Operator:** Marketing  
**Injection date:** 2005-10-14 15:34:53+07:00  
**Location:** Vial 4  
**Type:** Sample  
**Calib Level:**  
**Sample amount:** 0.00

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$$\text{Absolute Amt of } x = \text{Response}_x * RF_x * M * D$$

*Response<sub>x</sub>* is the response of peak X

*RF<sub>x</sub>* is the response factor for component X, (= Amount<sub>x</sub> / Response<sub>x</sub>)

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*Response<sub>x</sub>* is the area (or height) of peak X

*RF<sub>x</sub>* is the response factor

*M* is the multiplier

*D* is the dilution factor

$\sum(\text{Response} * RF)$  is the total of all the (RESPONSE)(RF) products for all peaks including peak x

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*Response Ratio* is  $\text{Response}_x / \text{Response}_{\text{ISTD}}$

*RF<sub>x</sub>* is the response factor for component X

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

Properties CDS Settings

Name:

Project folder path:  Select your project path

Include project groups in project path


Description:

Applications:  OpenLAB CDS  Don't unmark

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Edit System Settings

## Additional parameters:

You have chosen **Internal** as the authentication provider.  
Please specify the following additional parameters.

### Authentication Parameters

Create an account that will be used to administer the system.

#### Create Administrator Account

Name:	<input type="text" value="Create User Name"/>
Full Name:	<input type="text" value="Name for Reporting"/>
Password:	<input type="text" value="Create your Password"/>
Confirm password:	<input type="text" value="Confirm your password"/>

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
1. Double Click OpenLab Control Panel

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Name (ID): SP

Description:

General    Group Membership    **Role Membership**

Member of:

	Name	Description
<input checked="" type="checkbox"/>	Everything	All privileges
<input type="checkbox"/>	System Administrator	Manage users and security settings
<input type="checkbox"/>	Instrument Administrator	Manage instruments and locations
<input type="checkbox"/>	Project Administrator	Manage projects and project groups
<input type="checkbox"/>	Instrument User	View and run instruments
<input type="checkbox"/>	Technician	Laboratory technician
<input type="checkbox"/>	Chemist	Analytical chemist
<input type="checkbox"/>	admin_project	
<input type="checkbox"/>	admin_instrument	

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